



Product Overview

PRODUCT:	Recovered Human Plasma, Pooled
PART NUMBER:	5929
FORM / APPEARANCE:	Liquid, amber to straw, no grossly hemolyzed, lipemic, or icteric material
FILTRATION:	Filtration available
ADDITIVES:	None
STORAGE:	≤ -15°C
EXPIRATION:	5 year from date of manufacture

Analytical Properties

ANALYTE	TARGET
TOTAL PROTEIN:	5.7 g/dL
ALBUMIN:	3.4 g/dL
UREA:	13 mg/dL
CREATININE:	0.8 mg/dL
CALCIUM:	7.0 mg/dL
ALP:	49 U/L
AST:	15 U/L
POTASSIUM:	3.3 mmol/L
SODIUM:	167 mmol/L
CHLORIDE:	72 mmol/L
ALT:	17 U/L
TOTAL BILIRUBIN:	0.2 mg/dL
GLUCOSE:	389 mg/dL

FOR IN VITRO DIAGNOSTIC MANUFACTURING USE ONLY

All human plasma units used to prepare this product have been tested by FDA approved tests and found to be non-reactive for HbsAg, Anti-HIV I/II, Anti-HCV, HIV-1 RNA, and HCV RNA. Donors of human plasma units used in the making of this product were tested and found negative for Syphilis according to FDA guidelines.

CAUTION: FOR USE IN RESEARCH OR REMANUFACTURE OF NON-INJECTABLE PRODUCTS NOT SUBJECT TO LICENSE UNDER PHS ACT SECTION 351, ONLY.

CAUTION: ALL BLOOD PRODUCTS SHOULD BE HANDLED AND TREATED AS POTENTIALLY INFECTIOUS.