



Product Overview

PRODUCT:	Drug Free Urine, Unprocessed
PART NUMBER:	6907
FORM / APPEARANCE:	Liquid, light yellow
FILTRATION:	None
ADDITIVES:	None
STORAGE:	< - 20 °C
EXPIRATION:	3 years from date of manufacture

Analytical Properties

ANALYTE	VALUE
SPECIFIC GRAVITY:	≥ 1.010
GLUCOSE:	Negative
NITRITE:	Negative
PROTEIN:	Negative
KETONE:	Negative
BILIRUBIN:	Negative
UROBILINOGEN:	0.2 – 1.0 EU/dL
LEUKOCYTE:	Negative
BLOOD:	Negative
hCG:	As Found
ACETAMINOPHEN:	< 0.5 µg/mL
AMPHETAMINE:	< 50 ng/mL
BARBITURATES:	< 40 ng/mL
BENZODIAZEPINES:	Negative
COCAINE:	< 25 ng/mL
IBUPROFEN:	< 2 µg/mL
MARIJUANA:	< 2 ng/mL
METHADONE:	< 25 ng/mL
METHAMPHETAMINE:	< 50 ng/mL
METHAQUALONE:	< 25 ng/mL
METHYLENedioxy METHAMPHETAMINE:	< 40 ng/mL
6-MAM:	< 4 ng/mL
OPIATE:	Negative
OXYCODONE:	Negative
PCP:	< 2.5 ng/mL
PHENCYCLIDINE:	Negative
SALYCILATES:	< 5 µg/mL

TRICYCLIC ANTIDEPRESSANTS:	Negative
PROPOXYPHENE	< 75 ng/mL
NICOTINE:	0.0
ETHANOL:	< 5.0 mg/dL

FOR IN VITRO DIAGNOSTIC MANUFACTURING USE ONLY

Donors have been tested by FDA approved tests and found to be negative or non-reactive for HBsAg, Anti-HIV I/II, Anti-HCV, HIV-1 RNA, HCV RNA. Donors were tested and found negative for RPR according to FDA guidelines.

CAUTION: FOR LABORATORY REAGENT USE OR MANUFACTURING INTO NONINJECTABLE PRODUCTS ONLY. THE USE OF THIS MATERIAL IN HUMAN MEDICATION IS PROHIBITED.

