



Product Overview

| | |
|--------------------|-------------------------------------------------------------------------------|
| PRODUCT: | aHAV Negative, Delipidized Human Serum |
| PART NUMBER: | 5080 |
| FORM / APPEARANCE: | Liquid, amber to straw, free of visible foreign matter, clots or precipitates |
| FILTRATION: | 0.22 micron-filtered |
| ADDITIVES: | No preservatives |
| STORAGE: | < - 10 °C |
| EXPIRATION: | 5 years from date of manufacture |

Analytical Properties

| ANALYTE | VALUE |
|---------------------------|------------------|
| TOTAL PROTEIN: | 5.5 – 6.0 g/dL |
| aHAV: | ≤ 1.0 S/CO |
| TRIGLYCERIDE: | < 10 mg/dL |
| CHOLESTEROL: | < 10 mg/dL |
| CALCIUM: | < 10 mg/dL |
| SODIUM: | 130 – 148 mmol/L |
| CHLORIDE: | 90 – 110 mmol/L |
| POTASSIUM: | 3.2 – 5.4 mmol/L |
| AMMONIA: | < 1.0 mg/dL |
| pH: | 7.1 – 7.4 |
| TURBIDITY (O.D. @ 546nm): | < 0.14 |
| BIOBURDEN: | < 100 CFUs/mL |

FOR IN VITRO DIAGNOSTIC MANUFACTURING USE ONLY

All human plasma units used to prepare this product have been tested by FDA approved tests and found to be non-reactive for HbsAg, Anti-HIV I/II, Anti-HCV, HIV-1 RNA, and HCV RNA. Donors of human plasma units used in the making of this product were tested and found negative for Syphilis according to FDA guidelines.

CAUTION: FOR USE IN RESEARCH OR REMANUFACTURE OF NON-INJECTABLE PRODUCTS NOT SUBJECT TO LICENSE UNDER PHS ACT SECTION 351, ONLY.

CAUTION: ALL BLOOD PRODUCTS SHOULD BE HANDLED AND TREATED AS POTENTIALLY INFECTIOUS